

## ALS Airway CPR (AW101) with an Advanced Airway

This manikin has a high fidelity advanced airway providing feedback on how well you intubate. This improves cooperation when performing CPR with an advanced airway and while intubating on a moving chest. This product is also suitable for training in BLS skills.

### CPR with Advanced Airway

- Perform 30:2 BLS CPR
- Perform 10:1 Asynchronous CPR  
(Continuous Compressions With  
A Ventilation Every 6 Seconds)
- Visually See Interaction of  
Compressions and Ventilation



### Intubate

- Reports How Quickly Airway Device Established
- Reports Whether Cuff is Holding
- Reports Efficiency and Gives A Score
- Also Practice with Compressions Ongoing

### Team Work

- Feedback on Skills Working Together
- Score For the Total Performance
- Record Team or Individuals
- Recall Data For Team Debriefing

Note: Available in full body

★ Leading the World in Feedback that Creates Competence