

Baby CPR Pro (B201)

Advanced CPR Training

A full body baby manikin with an excellent look and feel. Accurate for even very low volumes of air. It features the award winning real time feedback system.

Real Time Feedback

Compressions - Depth, Rate, Release, Too Deep, Hitting Bottom

Ventilations - Volume and Rate

CPR - Suitable for two hand method and in arm method



See Chest Rise



Correct Head Position

Adjust Lung Capacity

Results Management

- Review and Print
- Class List with Scores
- Export Results To Excel

Focus On Ventilations

This product includes several special activities to improve ventilations on a baby. Respiratory difficulties occur more often with babies than with adults. These practice activities will improve how you use the BVM.

Rate & Volume Control

- Visual Display of Volume
- Running Timer for Air Flow

Lung Sizes

- Selection of 5 lungs sizes
- Single click to set size



Improve Tidal Flow

- Shows in and out movement of air
- Shows Proper Release of the BVM

Note: Manikin simulates 1 month to 12 month old baby

★ Leading the World in Feedback that Creates Competence