



one
size

TOPLINE PRODUCTS PROVIDE COMPRESSION,
WARMTH AND ENHANCED PROPRIOCEPTION.

- ORTHESES IN 2MM-THICK ASTRAKHAN NEOPRENE.
- DESIGNED FOR VARYING COMPRESSION LEVELS
ACCORDING TO USER'S NEEDS.
- LOW COMPRESSION FOR EVERYDAY USE.
- PAIN RELIEF.



ALL YOU NEED FOR PROTECTION IN SPORTS: You are likely to find all you need in our range of products: braces for the knee, wrist, elbow, thigh, calf, and back support, designed especially to provide specific solutions to each type of injury or pain.

THE SATISFACTION OF THINGS WELL DONE: Products designed and made in our SPAIN-based manufacturing facility, in compliance with strict quality control standards, which enables us to affix the CE marking to our neoprene items.

LATEX-FREE: Committed to prevent health problems related with exposure to latex, our full range of products are latex-free.

FOCUSED ON EVERY DETAIL: Distinctive features that will decide you to choose TopLine as your favourite brand:

- Strong fastening for perfect fit
- Inner towelling lining
- Comfortable edges
- 100% medical grade silicone inserts




Prim
Customers Service - Export
Pol. Ind. Nº1 c/C nº20. 28938 Móstoles - Madrid - SPAIN
Tel.: +34 91 334 25 32 - Fax: +34 91 334 25 65
export@ortoprim.com



one
size

Protection in Sports

NEOPRENE



BRACES FOR THE KNEE,
ANKLE, WRIST,
ELBOW, THIGH,
CALF AND BACK

latex
100% free



one
size



Wrap-around Wrist Support Brace

Easy-to-use, effective wrist support with adjustable strap for fast and easy fastening.

KEY FEATURES: Provides warmth, adjustable levels of compression and muscle tension release.

INDICATIONS: Arthrosis, arthritis, tendinitis, bursitis, mild joint traumas.

TLS265 ONE SIZE



Wrist Strap Bandage

With the same features as the TLS265, it restricts movement to relieve pain or feeling of discomfort in the CMC joint.

KEY FEATURES: Provides warmth, adjustable levels of compression and muscle tension release.

INDICATIONS: Basal thumb arthrosis, arthritis, tendinitis, bursitis, strains and mild trauma.

TLS264 ONE SIZE



Elbow Brace

Orthesis with two straps – one above and the other below the elbow – that ensure a very good fit. Specifically designed to allow normal elbow movement.

KEY FEATURES: Provides warmth, compression and muscle tension release.

INDICATIONS: Epicondylitis (tennis elbow), epitrochleitis (golf elbow) and post trauma.

TLS273 ONE SIZE



Epicondylitis Strap

Unsurpassed design for a band with a 100% medical grade silicone pad for perfect tension release and to prevent injuries caused by overexertion. The outer cinch strap comes with a hook-and-loop closure for secure fit. 100% silicone pad.

KEY FEATURES: Provides warmth, adjustable levels of compression and muscle tension release.

INDICATIONS: Epicondylitis (tennis elbow) and epitrochleitis (golf elbow).

TLS272 ONE SIZE



Thigh Support

One-size cone-shaped orthesis with two adjustable straps for variable levels of compression.

KEY FEATURES: Provides warmth, compression and muscle tension release.

INDICATIONS: Fibrillar fractures, muscle strain, post trauma.

TLS247 ONE SIZE



Open Patella Knee Brace

Orthesis with open knee and double layer to maintain the proper alignment of the patella. The spiral stays on each side of the brace provide additional support.

KEY FEATURES: Provides warmth, compression, and medial and lateral stability while keeping the proper alignment of the patella.

INDICATIONS: Joint pain, bursitis, tendinitis, arthritis, arthrosis, patella injuries, collateral ligament strains, Osgood-Schlatter disease, patellar chondromalacia, post trauma.

TLS218 ONE SIZE



Patella Strap

Simple, effective product that allows freedom of movement and has an adjustable compression zone for a perfect fit.

KEY FEATURES: Applies pressure on the tendon below the knee cap and mild compression at the patellar tendon.

INDICATIONS: Osgood-Schlatter disease, patellar chondromalacia, patellar tendinitis, sport that involves excessive knee extension.

TLS232 ONE SIZE



Calf Support

One-size orthesis with three adjustable straps for variable levels of compression.

KEY FEATURES: Provides warmth, compression and muscle tension release.

INDICATIONS: Fibrillar fractures, muscle strain, post trauma.

TLS220 ONE SIZE



Back Support Belt

Designed for everyday use when slight compression is required for unrestricted movement. 24cm height.

KEY FEATURES: Provides warmth, variable levels of compression and muscle tension release.

INDICATIONS: Mild low back pain, conditions that require light compression.

TLS255 ONE SIZE



Ankle Support Brace

Open heel for wide joint movement. The area of the brace in contact with the plantar surface is made with breathable material.

KEY FEATURES: Provides warmth, compression, medial and lateral stability and muscle tension release.

INDICATIONS: Mild strains and strain prevention during sports, tendinitis.

TLS212 ONE SIZE



Lumbosacral Back Support

Orthesis construction that provides moderate support. Reinforced by four stays. 24cm back height that reduces to 16cm at the front for easy fit.

KEY FEATURES: Provides warmth, variable levels of compression and muscle tension release.

INDICATIONS: Acute or chronic low back pain, minor trauma, low back pain during sports.

TLS256 ONE SIZE