


# TENS-100

## Single Channel Electrical Stimulator

### Features


#### ■ Six preprogrammed modes:

Tapping (slow), Tapping (fast), Kneading (slow), Kneading (fast), Pressing and Pressing (Continuous) modes are available.



#### ■ Comfortable new spike waveform:

The gentle, electrified spike waveform reduces unpleasant tingling and provides a pleasant, comfortable massage.



#### ■ Automatic timer:

The convenient automatic shutdown function turns the unit off after 15 minutes.

#### ■ Easy operation:

The large dial makes it easy to adjust intensity.

### Specifications

**Power source:** DC 4.5V, three AA dry-cell batteries

**Output current:** 10mA (rms) or less at 1k ohm load

**Output voltage:** 130V or less


**Pulse shape:** spike waveform

**Pulse rate:** 160Hz or less

**Dimensions:** 140(H) × 68(W) × 22(D) mm

**Weight:** 130g ± 20% including batteries

### Standard Accessories



**ITO** CO., LTD.

3-3-3 Toyotama-Minami, Nerima-ku, Tokyo 176-0014, Japan

TEL: 81-3-3994-4619 FAX: 81-3-3994-1465

URL: <http://www.itolator.co.jp/english/> E-Mail: [itocoltd@itolator.co.jp](mailto:itocoltd@itolator.co.jp)



ISO14001 (Toyotama factory)

\*Designs and specifications are subject to change without notice.

Printed in Japan No.178 D0701VF-1